



N4BP1 blocking peptide (CDBP5781)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	N4BP1 NEDD4 binding protein 1 [Homo sapiens (human)]
Official Symbol	N4BP1
Synonyms	N4BP1; NEDD4 binding protein 1; NEDD4-binding protein 1
Entrez Gene ID	9683
mRNA Refseq	NM_153029
Protein Refseq	NP_694574
UniProt ID	O75113
Function	protein binding