



NGFRAP1 blocking peptide (CDBP5782)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NGFRAP1 nerve growth factor receptor (TNFRSF16) associated protein 1 [Homo sapiens (human)]
Official Symbol	NGFRAP1
Synonyms	NGFRAP1; nerve growth factor receptor (TNFRSF16) associated protein 1; Bex; BEX3; NADE; HGR74; DXS6984E; protein BEX3; brain-expressed X-linked protein 3; p75NTR-associated cell death executor; ovarian granulosa cell 13.0 kDa protein HGR74
Entrez Gene ID	27018
mRNA Refseq	NM_001282674
Protein Refseq	NP_001269603
UniProt ID	Q00994

Pathway	Cell death signalling via NRAGE; NADE modulates death signalling; Neurotrophin signaling pathway; Signal Transduction; Signalling by NGF; p75 NTR receptor-mediated signalling; p75(NTR)-mediated signaling
Function	cysteine-type endopeptidase activator activity involved in apoptotic process; death receptor binding; metal ion binding; protein binding