



NANOS1 blocking peptide (CDBP5798)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NANOS1 nanos homolog 1 (Drosophila) [Homo sapiens (human)]
Official Symbol	NANOS1
Synonyms	NANOS1; nanos homolog 1 (Drosophila); NOS1; SPGF12; nanos homolog 1; NOS-1; EC_Rep1a
Entrez Gene ID	340719
mRNA Refseq	NM_199461
Protein Refseq	NP_955631
UniProt ID	Q8WY41
Function	RNA binding; protein binding; translation repressor activity; zinc ion binding