



NANOS2 blocking peptide (CDBP5801)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NANOS2 nanos homolog 2 (Drosophila) [Homo sapiens (human)]
Official Symbol	NANOS2
Synonyms	NANOS2; nanos homolog 2 (Drosophila); NOS2; nanos homolog 2; NOS-2
Entrez Gene ID	339345
mRNA Refseq	NM_001029861
Protein Refseq	NP_001025032
UniProt ID	P60321
Function	RNA binding; zinc ion binding