



NANOS3 blocking peptide (CDBP5802)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NANOS3 nanos homolog 3 (Drosophila) [Homo sapiens (human)]
Official Symbol	NANOS3
Synonyms	NANOS3; nanos homolog 3 (Drosophila); NOS3; NANOS1L; nanos homolog 3
Entrez Gene ID	342977
mRNA Refseq	NM_001098622
Protein Refseq	NP_001092092
UniProt ID	P60323
Function	RNA binding; zinc ion binding