



NAA40 blocking peptide (CDBP5803)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NAA40 N(alpha)-acetyltransferase 40, NatD catalytic subunit [Homo sapiens (human)]
Official Symbol	NAA40
Synonyms	NAA40; N(alpha)-acetyltransferase 40, NatD catalytic subunit; NAT11; PATT1; N-alpha-acetyltransferase 40; natD catalytic subunit; N-alpha acetyl transferase 40; N-acetyltransferase 11 (GCN5-related, putative); N-alpha-acetyltransferase 40, NatD catalytic subunit; N(alpha)-acetyltransferase 40, NatD catalytic subunit, homolog
Entrez Gene ID	79829
mRNA Refseq	NM_001300800
Protein Refseq	NP_001287729
UniProt ID	Q86UY6

Pathway	metapathway biotransformation
Function	N-acetyltransferase activity