



NCLN blocking peptide (CDBP5804)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NCLN nicalin [Homo sapiens (human)]
Official Symbol	NCLN
Synonyms	NCLN; nicalin; NET59; nicastrin-like protein; nicalin homolog (zebrefish)
Entrez Gene ID	56926
mRNA Refseq	NM_020170
Protein Refseq	NP_064555
UniProt ID	Q969V3
Function	peptidase activity; protein binding