



GBAS blocking peptide (CDBP5820)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a member of the NipSnap family of proteins that may be involved in vesicular transport. The encoded protein is localized to mitochondria and plays a role in oxidative phosphorylation. A pseudogene of this gene is located on the long arm of chromosome 2. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Feb 2011]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	GBAS glioblastoma amplified sequence [Homo sapiens (human)]
Official Symbol	GBAS
Synonyms	GBAS; glioblastoma amplified sequence; NIPSNAP2; protein NipSnap homolog 2; 4-nitrophenylphosphatase domain and non-neuronal SNAP25-like 2
Entrez Gene ID	2631
mRNA Refseq	NM_001202469

Protein Refseq	NP_001189398
UniProt ID	O75323
Function	protein binding