



NIPSNAP3A blocking peptide (CDBP5822)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NIPSNAP3A belongs to a family of proteins with putative roles in vesicular transport (Buechler et al., 2004 [PubMed 15177564]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NIPSNAP3A nipsnap homolog 3A (C. elegans) [Homo sapiens (human)]
Official Symbol	NIPSNAP3A
Synonyms	NIPSNAP3A; nipsnap homolog 3A (C. elegans); TASSC; HSPC299; NIPSNAP4; protein NipSnap homolog 3A; protein NipSnap homolog 4; target for Salmonella secreted protein C
Entrez Gene ID	25934
mRNA Refseq	NM_015469
Protein Refseq	NP_056284

UniProt ID	Q9UFN0
Function	protein binding