



NIPSNAP3B blocking peptide (CDBP5823)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NIPSNAP3B belongs to a family of proteins with putative roles in vesicular trafficking (Buechler et al., 2004 [PubMed 15177564]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NIPSNAP3B nipsnap homolog 3B (C. elegans) [Homo sapiens (human)]
Official Symbol	NIPSNAP3B
Synonyms	NIPSNAP3B; nipsnap homolog 3B (C. elegans); NIPSNAP3; protein NipSnap homolog 3B; SNAP1
Entrez Gene ID	55335
mRNA Refseq	NM_018376
Protein Refseq	NP_060846

