



NKX2-4 blocking peptide (CDBP5828)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NKX2-4 NK2 homeobox 4 [Homo sapiens (human)]
Official Symbol	NKX2-4
Synonyms	NKX2-4; NK2 homeobox 4; NKX2D; NKX2.4; homeobox protein Nkx-2.4; NK2 transcription factor homolog D; NK2 transcription factor related, locus 4
Entrez Gene ID	644524
mRNA Refseq	NM_033176
Protein Refseq	NP_149416
UniProt ID	Q9H2Z4
Function	molecular_function; sequence-specific DNA binding; sequence-specific DNA binding

