



NKX2-8 blocking peptide (CDBP5830)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NKX2-8 NK2 homeobox 8 [Homo sapiens (human)]
Official Symbol	NKX2-8
Synonyms	NKX2-8; NK2 homeobox 8; NKX2H; NKX2.8; Nkx2-9; homeobox protein Nkx-2.8; NK-2 homolog 8; NK-2 homolog H; homeobox protein NK-2 homolog H; NK2 transcription factor related, locus 8
Entrez Gene ID	26257
mRNA Refseq	NM_014360
Protein Refseq	NP_055175
UniProt ID	O15522

Function	double-stranded DNA binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity
-----------------	--
