



Mrub_1588 blocking peptide (CDBP5831)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Mrub_1588 hypothetical protein [Meiothermus ruber DSM 1279]
Official Symbol	Mrub_1588
Synonyms	Mrub_1588; hypothetical protein
Entrez Gene ID	8879693