



NOL3 blocking peptide (CDBP5840)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes an anti-apoptotic protein that has been shown to down-regulate the enzyme activities of caspase 2, caspase 8 and tumor protein p53. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NOL3 nucleolar protein 3 (apoptosis repressor with CARD domain) [Homo sapiens (human)]
Official Symbol	NOL3
Synonyms	NOL3; nucleolar protein 3 (apoptosis repressor with CARD domain); ARC; FCM; MYP; NOP; NOP30; nucleolar protein 3; nucleolar protein of 30 kDa; muscle-enriched cytoplasmic protein
Entrez Gene ID	8996
mRNA Refseq	NM_001185057
Protein Refseq	NP_001171986

UniProt ID	O60936
Function	RNA binding; identical protein binding; protein binding