



Pmaip1 blocking peptide (CDBP5843)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Pmaip1 phorbol-12-myristate-13-acetate-induced protein 1 [Mus musculus (house mouse)]
Official Symbol	Pmaip1
Synonyms	Pmaip1; phorbol-12-myristate-13-acetate-induced protein 1; Noxa; protein Noxa
Entrez Gene ID	58801
mRNA Refseq	NM_021451
Protein Refseq	NP_067426
UniProt ID	Q9JM54
Pathway	Activation of BH3-only proteins; Activation of NOXA and translocation to mitochondria; Apoptosis; BH3-only proteins associate with and inactivate anti-apoptotic BCL-2 members;

Intrinsic Pathway for Apoptosis; Viral carcinogenesis; miRNA regulation of DNA Damage Response; p53 signaling pathway

Function	protein binding
-----------------	-----------------
