



WDR74 blocking peptide (CDBP5850)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	WDR74 WD repeat domain 74 [Homo sapiens (human)]
Official Symbol	WDR74
Synonyms	WDR74; WD repeat domain 74; WD repeat-containing protein 74; NOP seven-associated protein 1
Entrez Gene ID	54663
mRNA Refseq	NM_018093
Protein Refseq	NP_060563
UniProt ID	Q6RFH5
Function	molecular_function