



NSRP1 blocking peptide (CDBP5851)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NSRP1 nuclear speckle splicing regulatory protein 1 [Homo sapiens (human)]
Official Symbol	NSRP1
Synonyms	NSRP1; nuclear speckle splicing regulatory protein 1; CCDC55; NSrp70; HSPC095; coiled-coil domain containing 55; nuclear speckle-related protein 70; coiled-coil domain-containing protein 55
Entrez Gene ID	84081
mRNA Refseq	NM_001261467
Protein Refseq	NP_001248396
UniProt ID	Q9H0G5

Function

mRNA binding; poly(A) RNA binding; protein binding
