



OCIAD1 blocking peptide (CDBP5859)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	OCIAD1 OCIA domain containing 1 [Homo sapiens (human)]
Official Symbol	OCIAD1
Synonyms	OCIAD1; OCIA domain containing 1; OCIA; ASRIJ; TPA018; OCIA domain-containing protein 1; ovarian carcinoma immunoreactive antigen
Entrez Gene ID	54940
mRNA Refseq	NM_001079839
Protein Refseq	NP_001073308
UniProt ID	Q9NX40
Function	protein binding