



ORAI2 blocking peptide (CDBP5869)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ORAI2 ORAI calcium release-activated calcium modulator 2 [Homo sapiens (human)]
Official Symbol	ORAI2
Synonyms	ORAI2; ORAI calcium release-activated calcium modulator 2; CBCIP2; C7orf19; MEM142B; TMEM142B; protein orai-2; H_NH0514P08.8; putative protein ORAI2-2; transmembrane protein 142B; CAP-binding protein complex interacting protein 2; CAP-binding protein complex-interacting protein 2
Entrez Gene ID	80228
mRNA Refseq	NM_001126340
Protein Refseq	NP_001119812
UniProt ID	Q96SN7

Pathway	Calcium signaling pathway
Function	protein binding