



ORAI3 blocking peptide (CDBP5870)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ORAI3 ORAI calcium release-activated calcium modulator 3 [Homo sapiens (human)]
Official Symbol	ORAI3
Synonyms	ORAI3; ORAI calcium release-activated calcium modulator 3; TMEM142C; protein orai-3; transmembrane protein 142C
Entrez Gene ID	93129
mRNA Refseq	NM_152288
Protein Refseq	NP_689501
UniProt ID	Q9BRQ5
Pathway	Calcium signaling pathway

Functionprotein binding
