



ORMDL1 blocking peptide (CDBP5872)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ORMDL1 ORMDL sphingolipid biosynthesis regulator 1 [Homo sapiens (human)]
Official Symbol	ORMDL1
Synonyms	ORMDL1; ORMDL sphingolipid biosynthesis regulator 1; ORM1-like protein 1; adoplín-1
Entrez Gene ID	94101
mRNA Refseq	NM_001128150
Protein Refseq	NP_001121622
UniProt ID	Q9P0S3