



OTUD4 blocking peptide (CDBP5873)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Alternatively spliced transcript variants have been found for this gene. The smaller protein isoform encoded by the shorter transcript variant is found only in HIV-1 infected cells. [provided by RefSeq, Jul 2010]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	OTUD4 OTU deubiquitinase 4 [Homo sapiens (human)]
Official Symbol	OTUD4
Synonyms	OTUD4; OTU deubiquitinase 4; HIN1; DUBA6; HSHIN1; OTU domain-containing protein 4; OTU domain containing 4; HIV-1 induced protein HIN-1; HIV-1-induced protein HIN-1
Entrez Gene ID	54726
mRNA Refseq	NM_001102653
Protein Refseq	NP_001096123

UniProt ID	Q01804
Function	poly(A) RNA binding; ubiquitin-specific protease activity