



TP53INP1 blocking peptide (CDBP5877)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TP53INP1 tumor protein p53 inducible nuclear protein 1 [Homo sapiens (human)]
Official Symbol	TP53INP1
Synonyms	TP53INP1; tumor protein p53 inducible nuclear protein 1; SIP; Teap; p53DINP1; TP53DINP1; TP53INP1A; TP53INP1B; tumor protein p53-inducible nuclear protein 1; p53-inducible p53DINP1; stress-induced protein; p53-dependent damage-inducible nuclear protein 1
Entrez Gene ID	94241
mRNA Refseq	NM_001135733
Protein Refseq	NP_001129205
UniProt ID	Q96A56

Pathway	Direct p53 effectors; HTLV-I infection
Function	antioxidant activity; protein binding