



PALM3 blocking peptide (CDBP5885)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PALM3 paralemmin 3 [Homo sapiens (human)]
Official Symbol	PALM3
Synonyms	PALM3; paralemmin 3; paralemmin-3
Entrez Gene ID	342979
mRNA Refseq	NM_001145028
Protein Refseq	NP_001138500
UniProt ID	A6NDB9
Function	ATP binding; protein binding