



CUL9 blocking peptide (CDBP5887)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CUL9 cullin 9 [Homo sapiens (human)]
Official Symbol	CUL9
Synonyms	CUL9; cullin 9; PARC; H7AP1; cullin-9; CUL-9; UbcH7-associated protein 1; parkin-like cytoplasmic p53 binding protein; p53-associated parkin-like cytoplasmic protein
Entrez Gene ID	23113
mRNA Refseq	NM_015089
Protein Refseq	NP_055904
UniProt ID	Q8IWT3
Function	ATP binding; protein binding; ubiquitin protein ligase binding; zinc ion binding