



PDCD1LG2 blocking peptide (CDBP5896)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PDCD1LG2 programmed cell death 1 ligand 2 [Homo sapiens (human)]
Official Symbol	PDCD1LG2
Synonyms	PDCD1LG2; programmed cell death 1 ligand 2; B7DC; Btdc; PDL2; CD273; PD-L2; PDCD1L2; bA574F11.2; B7-DC; PD-1 ligand 2; PD-1-ligand 2; PDCD1 ligand 2; butyrophilin B7-DC; programmed death ligand 2; B7 dendritic cell molecule
Entrez Gene ID	80380
mRNA Refseq	NM_025239
Protein Refseq	NP_079515
UniProt ID	Q9BQ51

Pathway	Adaptive Immune System; Cell adhesion molecules (CAMs); Costimulation by the CD28 family; Immune System; PD-1 signaling
Function	molecular_function