



PEAR1 blocking peptide (CDBP5897)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PEAR1 is a platelet receptor that signals upon the formation of platelet-platelet contacts independent of platelet activation and secondary to platelet aggregation (Nanda et al., 2005 [PubMed 15851471]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PEAR1 platelet endothelial aggregation receptor 1 [Homo sapiens (human)]
Official Symbol	PEAR1
Synonyms	PEAR1; platelet endothelial aggregation receptor 1; JEDI; MEGF12; hPEAR1; multiple EGF-like-domains 12; multiple EGF-like domains protein 12; multiple epidermal growth factor-like domains protein 12
Entrez Gene ID	375033
mRNA Refseq	NM_001080471

Protein Refseq [NP_001073940](#)

UniProt ID Q5VY43
