



PELI1 blocking peptide (CDBP5898)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PELI1 pellino E3 ubiquitin protein ligase 1 [Homo sapiens (human)]
Official Symbol	PELI1
Synonyms	PELI1; pellino E3 ubiquitin protein ligase 1; E3 ubiquitin-protein ligase pellino homolog 1; pellino-1; pellino homolog 1; protein pellino homolog 1; pellino-related intracellular-signaling molecule
Entrez Gene ID	57162
mRNA Refseq	NM_020651
Protein Refseq	NP_065702
UniProt ID	Q96FA3

Pathway	Activated TLR4 signalling; Cytokine Signaling in Immune system; IRAK1 recruits IKK complex; IRAK1 recruits IKK complex upon TLR7/8 or 9 stimulation; Immune System; Innate Immune System; Interleukin-1 signaling; MyD88 cascade initiated on plasma membrane
Function	ligase activity; protein binding; ubiquitin-protein transferase activity