



FAM120B blocking peptide (CDBP5903)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FAM120B family with sequence similarity 120B [Homo sapiens (human)]
Official Symbol	FAM120B
Synonyms	FAM120B; family with sequence similarity 120B; CCPG; PGCC1; KIAA1838; dJ894D12.1; constitutive coactivator of peroxisome proliferator-activated receptor gamma; PPARG constitutive coactivator 1; PPARGgamma constitutive coactivator 1; constitutive coactivator of PPAR-gamma
Entrez Gene ID	84498
mRNA Refseq	NM_001286379
Protein Refseq	NP_001273308
UniProt ID	Q96EK7

