



ANP32E blocking peptide (CDBP5905)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ANP32E acidic (leucine-rich) nuclear phosphoprotein 32 family, member E [Homo sapiens (human)]
Official Symbol	ANP32E
Synonyms	ANP32E; acidic (leucine-rich) nuclear phosphoprotein 32 family, member E; LANPL; LANP-L; acidic leucine-rich nuclear phosphoprotein 32 family member E; LANP-like protein; leucine-rich acidic nuclear protein like
Entrez Gene ID	81611
mRNA Refseq	NM_001136478
Protein Refseq	NP_001129950
UniProt ID	Q9BTT0

Function

histone binding; phosphatase inhibitor activity
