



PHC3 blocking peptide (CDBP5906)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PHC3 polyhomeotic homolog 3 (Drosophila) [Homo sapiens (human)]
Official Symbol	PHC3
Synonyms	PHC3; polyhomeotic homolog 3 (Drosophila); EDR3; HPH3; polyhomeotic-like protein 3; polyhomeotic like 3; homolog of polyhomeotic 3; early development regulator 3; early development regulatory protein 3
Entrez Gene ID	80012
mRNA Refseq	NM_024947
Protein Refseq	NP_079223
UniProt ID	Q8NDX5

Pathway	Cellular Senescence; Cellular responses to stress; Oxidative Stress Induced Senescence
Function	DNA binding; zinc ion binding