



PHF14 blocking peptide (CDBP5907)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PHF14 PHD finger protein 14 [Homo sapiens (human)]
Official Symbol	PHF14
Synonyms	PHF14; PHD finger protein 14
Entrez Gene ID	9678
mRNA Refseq	NM_014660
Protein Refseq	NP_055475
UniProt ID	O94880
Function	zinc ion binding