



PHLPP2 blocking peptide (CDBP5909)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PHLPP2 PH domain and leucine rich repeat protein phosphatase 2 [Homo sapiens (human)]
Official Symbol	PHLPP2
Synonyms	PHLPP2; PH domain and leucine rich repeat protein phosphatase 2; PHLPP1; PH domain leucine-rich repeat-containing protein phosphatase 2
Entrez Gene ID	23035
mRNA Refseq	NM_001289003
Protein Refseq	NP_001275932
UniProt ID	Q6ZVD8
Pathway	Adaptive Immune System; DAP12 interactions; DAP12 signaling; Disease; Downstream signal

transduction; Downstream signaling events of B Cell Receptor (BCR); Downstream signaling of activated FGFR; Fc epsilon receptor (FCERI) signaling

Function

metal ion binding; protein serine/threonine phosphatase activity
