



PIBF1 blocking peptide (CDBP5915)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PIBF1 progesterone immunomodulatory binding factor 1 [Homo sapiens (human)]
Official Symbol	PIBF1
Synonyms	PIBF1; progesterone immunomodulatory binding factor 1; PIBF; CEP90; C13orf24; progesterone-induced-blocking factor 1; progesterone-induced blocking factor 1
Entrez Gene ID	10464
mRNA Refseq	NM_006346
Protein Refseq	NP_006337
UniProt ID	Q8WXW3
Function	protein binding