



PLAC1 blocking peptide (CDBP5929)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PLAC1 placenta-specific 1 [Homo sapiens (human)]
Official Symbol	PLAC1
Synonyms	PLAC1; placenta-specific 1; CT92; OOSP2L; placenta-specific protein 1; cancer/testis antigen 92
Entrez Gene ID	10761
mRNA Refseq	NM_021796
Protein Refseq	NP_068568
UniProt ID	Q9HBJ0
Function	molecular_function