



PLAC4 blocking peptide (CDBP5932)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PLAC4 placenta-specific 4 [Homo sapiens (human)]
Official Symbol	PLAC4
Synonyms	PLAC4; placenta-specific 4; PRED78; C21orf115; placenta-specific protein 4
Entrez Gene ID	191585
mRNA Refseq	NM_182832
Protein Refseq	NP_878252
UniProt ID	Q8WY50