



PLAGL2 blocking peptide (CDBP5933)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Pleiomorphic adenoma gene-like 2 is a zinc-finger protein that recognizes DNA and/or RNA. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PLAGL2 pleiomorphic adenoma gene-like 2 [Homo sapiens (human)]
Official Symbol	PLAGL2
Synonyms	PLAGL2; pleiomorphic adenoma gene-like 2; ZNF900; zinc finger protein PLAGL2; C2H2-type zinc finger protein; pleiomorphic adenoma-like protein 2
Entrez Gene ID	5326
mRNA Refseq	NM_002657
Protein Refseq	NP_002648

UniProt ID	Q9UPG8
Function	RNA polymerase II core promoter proximal region sequence-specific DNA binding transcription factor activity involved in positive regulation of transcription; metal ion binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity