



PLEKHM2 blocking peptide (CDBP5935)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PLEKHM2 pleckstrin homology domain containing, family M (with RUN domain) member 2 [Homo sapiens (human)]
Official Symbol	PLEKHM2
Synonyms	PLEKHM2; pleckstrin homology domain containing, family M (with RUN domain) member 2; SKIP; pleckstrin homology domain-containing family M member 2; sifA and kinesin-interacting protein; PH domain-containing family M member 2; novel RUN and PH domain-containing protein; salmonella-induced filaments A and kinesin-interacting protein; SifA (Salmonella-induced filaments A) and kinesin-interacting protein
Entrez Gene ID	23207
mRNA Refseq	NM_015164
Protein Refseq	NP_055979

UniProt ID	Q8IWE5
Pathway	Salmonella infection
Function	kinesin binding; protein binding