



PLEKHM3 blocking peptide (CDBP5936)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PLEKHM3 pleckstrin homology domain containing, family M, member 3 [Homo sapiens (human)]
Official Symbol	PLEKHM3
Synonyms	PLEKHM3; pleckstrin homology domain containing, family M, member 3; DAPR; PLEKHM1L; pleckstrin homology domain-containing family M member 3; differentiation associated protein; PH domain-containing family M member 3
Entrez Gene ID	389072
mRNA Refseq	NM_001080475
Protein Refseq	NP_001073944
UniProt ID	Q6ZWE6

Functionmetal ion binding
