



## PLXDC2 blocking peptide (CDBP5937)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

### GENE INFORMATION

|                 |                                                                                                                                          |
|-----------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">PLXDC2 plexin domain containing 2 [ Homo sapiens (human) ]</a>                                                               |
| Official Symbol | PLXDC2                                                                                                                                   |
| Synonyms        | PLXDC2; plexin domain containing 2; TEM7R; plexin domain-containing protein 2; 1200007L24Rik; tumor endothelial marker 7-related protein |
| Entrez Gene ID  | <a href="#">84898</a>                                                                                                                    |
| mRNA Refseq     | <a href="#">NM_001282736</a>                                                                                                             |
| Protein Refseq  | <a href="#">NP_001269665</a>                                                                                                             |
| UniProt ID      | Q6UX71                                                                                                                                   |
| Function        | protein binding; receptor activity                                                                                                       |