



PLXDC2 blocking peptide (CDBP5938)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PLXDC2 plexin domain containing 2 [Homo sapiens (human)]
Official Symbol	PLXDC2
Synonyms	PLXDC2; plexin domain containing 2; TEM7R; plexin domain-containing protein 2; 1200007L24Rik; tumor endothelial marker 7-related protein
Entrez Gene ID	84898
mRNA Refseq	NM_001282736
Protein Refseq	NP_001269665
UniProt ID	Q6UX71
Function	protein binding; receptor activity