



DESI2 blocking peptide (CDBP5939)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	DESI2 desumoylating isopeptidase 2 [Homo sapiens (human)]
Official Symbol	DESI2
Synonyms	DESI2; desumoylating isopeptidase 2; DESI; DESI1; DeSI-2; PNAS-4; PPPDE1; CGI-146; FAM152A; C1orf121; desumoylating isopeptidase 1; PPPDE peptidase domain-containing protein 1; family with sequence similarity 152, member A
Entrez Gene ID	51029
mRNA Refseq	NM_001297746
Protein Refseq	NP_001284675
UniProt ID	Q9BSY9

Functionpeptidase activity
