



PPAPDC1A blocking peptide (CDBP5943)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PPAPDC1A phosphatidic acid phosphatase type 2 domain containing 1A [Homo sapiens (human)]
Official Symbol	PPAPDC1A
Synonyms	PPAPDC1A; phosphatidic acid phosphatase type 2 domain containing 1A; DPPL2; PPAPDC1; phosphatidate phosphatase PPAPDC1A; diacylglycerol pyrophosphate like 2; phosphatidic acid phosphatase type 2 domain-containing protein 1A
Entrez Gene ID	196051
mRNA Refseq	NM_001030059
Protein Refseq	NP_001025230
UniProt ID	Q5VZY2

Pathway	Fcγ receptor (FCGR) dependent phagocytosis; Immune System; Innate Immune System; Role of phospholipids in phagocytosis; triacylglycerol biosynthesis
Function	phosphatidate phosphatase activity