



PPAPDC1B blocking peptide (CDBP5944)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PPAPDC1B phosphatidic acid phosphatase type 2 domain containing 1B [Homo sapiens (human)]
Official Symbol	PPAPDC1B
Synonyms	PPAPDC1B; phosphatidic acid phosphatase type 2 domain containing 1B; DPPL1; HTPAP; phosphatidate phosphatase PPAPDC1B; diacylglycerol pyrophosphate like 1; phosphatidic acid phosphatase type 2 domain-containing protein 1B
Entrez Gene ID	84513
mRNA Refseq	NM_001102559
Protein Refseq	NP_001096029
UniProt ID	Q8NEB5

Pathway	Fc gamma receptor (FCGR) dependent phagocytosis; Immune System; Innate Immune System; Role of phospholipids in phagocytosis
Function	phosphatidate phosphatase activity