



PPAPDC2 blocking peptide (CDBP5945)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PPAPDC2 phosphatidic acid phosphatase type 2 domain containing 2 [Homo sapiens (human)]
Official Symbol	PPAPDC2
Synonyms	PPAPDC2; phosphatidic acid phosphatase type 2 domain containing 2; PDP1; PSDP; bA6J24.6; presqualene diphosphate phosphatase; PPAP2 domain-containing protein 2; polyisoprenoid diphosphate phosphatase type 1; phosphatidic acid phosphatase type 2 domain-containing protein 2
Entrez Gene ID	403313
mRNA Refseq	NM_203453
Protein Refseq	NP_982278

UniProt ID	Q8IY26
Function	hydrolase activity