



PPAPDC3 blocking peptide (CDBP5946)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PPAPDC3 phosphatidic acid phosphatase type 2 domain containing 3 [Homo sapiens (human)]
Official Symbol	PPAPDC3
Synonyms	PPAPDC3; phosphatidic acid phosphatase type 2 domain containing 3; NET39; C9orf67; probable lipid phosphate phosphatase PPAPDC3; phosphatidic acid phosphatase type 2 domain-containing protein 3
Entrez Gene ID	84814
mRNA Refseq	NM_032728
Protein Refseq	NP_116117
UniProt ID	Q8NBV4

Function	hydrolase activity
----------	--------------------
