



PRICKLE4 blocking peptide (CDBP5956)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C6ORF49 is a member of the LIM domain protein family (Teufel et al., 2005 [PubMed 15702247]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PRICKLE4 prickle homolog 4 (Drosophila) [Homo sapiens (human)]
Official Symbol	PRICKLE4
Synonyms	PRICKLE4; prickle homolog 4 (Drosophila); OBTP; OEBT; TOMM6; C6orf49; prickle-like protein 4; overexpressed breast tumor protein; over-expressed breast tumor protein; translocase of outer mitochondrial membrane 6 homolog
Entrez Gene ID	29964
mRNA Refseq	NM_013397
Protein Refseq	NP_037529

UniProt ID	Q2TBC4
Function	molecular_function; zinc ion binding