



LZTS3 blocking peptide (CDBP5957)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LZTS3 leucine zipper, putative tumor suppressor family member 3 [Homo sapiens (human)]
Official Symbol	LZTS3
Synonyms	LZTS3; PROSAPIP1; leucine zipper putative tumor suppressor 3; proSAP-interacting protein 1; ProSAP/Shank-interacting protein 1; proline rich synapse associated protein interacting protein 1
Entrez Gene ID	9762
mRNA Refseq	NM_001282533
Protein Refseq	NP_001269462
UniProt ID	O60299