



TBKBP1 blocking peptide (CDBP5958)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TBKBP1 is an adaptor protein that binds to TBK1 (MIM 604834) and is part of the interaction network in the TNF (MIM 191160)/NFKB (see MIM 164011) pathway (Bouwmeester et al., 2004 [PubMed 14743216]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TBKBP1 TBK1 binding protein 1 [Homo sapiens (human)]
Official Symbol	TBKBP1
Synonyms	TBKBP1; TBK1 binding protein 1; SINTBAD; ProSAPiP2; TANK-binding kinase 1-binding protein 1; TBK1-binding protein 1
Entrez Gene ID	9755
mRNA Refseq	NM_014726
Protein Refseq	NP_055541

UniProt ID	A7MCY6
Pathway	RIG-I-like receptor signaling pathway
Function	protein binding