



PRR5L blocking peptide (CDBP5959)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PRR5L proline rich 5 like [Homo sapiens (human)]
Official Symbol	PRR5L
Synonyms	PRR5L; proline rich 5 like; PROTOR2; proline-rich protein 5-like; protor-2; Protein observed with Rictor-2
Entrez Gene ID	79899
mRNA Refseq	NM_001160167
Protein Refseq	NP_001153639
UniProt ID	Q6MZQ0
Pathway	TOR signaling

Functionubiquitin protein ligase binding
